

http://www.setabiomedicals.com e-mail: info@setabiomedicals.com

Product number: K4-225 Product name: Seta-385-NHS

General Data

Molecular Mass:	528.4
Solubility:	Water, alcohol, DMF, DMSO
Insoluble:	Acetone, toluene
Storage:	Store in absence of light, desiccated and refrigerate

Description

• Hydrophilic, positively charged, amine-reactive label containing one NHS-ester group.

Applications

- Covalent labeling of proteins, amino-modified DNA and amino-modified oligonucleotides
- Resonance Energy Transfer (RET)
- Flow Cytometry
- Immunofluorescence
- Gene Expression
- Homogeneous Assays
- Assessment of protein structure

Advantages

- · Perfectly suited for excitation with the 330 and 380-nm diode lasers and UV light
- Extremely large Stokes' shift of over 120 nm
- High extinction coefficients and high quantum yields of 78%
- Low non-specific binding
- pH-insensitive between pH 3 and pH 10
- Good aqueous solubility; this label does not alter the solubility of the protein conjugate
- High photostability; e.g. compared to fluorescein
- Low molecular weight Seta dyes do not add substantial mass to the conjugates
- · Ideal for non-radioactive labeling of proteins, amino-modified DNA probes and amino-modified oligonucleotides

Spectral Data


Solvent System: phosphate buffer pH 7.4


Sample	Dye-to-protein Ratio	Absorption max. [nm]	Extinction Coefficient [M ⁻¹ cm ⁻¹]	Fluorescence max. [nm]	Q.Y. ¹ [%]
Free dye	—	384	25,000	508	78
IgG conjugate 1	2.0	387		505	44
IgG conjugate 2	4.0	387		505	35
IgG conjugate 3	6.0	387		505	29


¹Excitation at 380 nm




http://www.setabiomedicals.com e-mail: info@setabiomedicals.com







Product name: Seta-385-NHS



Total brightness (QY $\times \epsilon \times D/P$) vs. dye-to-protein ratio (D/P) of K4-225 - IgG conjugates in phosphate buffer (pH 7.4)