

Product number: K8-1377

Product name: Seta-690-NH-di-NHS

Date: 12/16/08 Page 1

General Data

Molecular Mass: 1259.28

Solubility: Water, Alcohol, DMF, DMSO **Insoluble:** Acetone, Chloroform, Toluene

Storage: Store out of light, desiccated and refrigerate

Description

Amine-reactive fluorescent label containing two reactive NHS-ester group.

Applications

- Covalent labeling of proteins, amino-modified DNA and amino-modified oligonucleotides
- Fluorescence Lifetime Label this label exhibits a distinct lifetime change upon binding to a biomolecule
- Resonance Energy Transfer (RET)
- Flow Cytometry
- Immunofluorescence
- · Gene Expression
- Homogeneous Assays
- · Assessment of protein structure

Advantages


- Perfectly suited for excitation with the 404-nm and 670-nm diode lasers, and UV light
- Sensitive; high extinction coefficients and high quantum yields up to 20% after covalent attachment to proteins
- Low non-specific binding
- $\bullet~$ pH-sensitive between pH 7 and pH 11 and pH-insensitive between pH 3 and pH 7
- Good aqueous solubility; this label does not alter the solubility of the protein conjugate
- High photostability; e.g. compared to fluorescein or Cy5TM
- Low molecular weight Seta dyes do not add substantial mass to the conjugates
- Ideal for non-radioactive labeling of proteins, amino-modified DNA probes and amino-modified oligonucleotides

Spectral Data

Solvent System: phosphate buffer pH 6.0

Sample	Dye-to-protein Ratio	Absorption max. [nm]	Extinction Coefficient [M ⁻¹ ·cm ⁻¹]	Fluorescence* max. [nm]	Quantum Yield [%]
Free dye	_	693	113,000	714	14
BSA conjugate 1	0.5	711		726	17
BSA conjugate 2	1.0	711		726	14
BSA conjugate 3	1.5	710		728	11
BSA conjugate 4	2.0	710		730	10


^{*} Excitation at 645 nm







Product number: K8-1377

Product name: Seta-690-NH-di-NHS







Absorption and emission spectrum of **K8-1377** — **BSA conjugate (D/P = 1.5)** in phosphate buffer (pH 6.0)

