



http://www.setabiomedicals.com e-mail: info@setabiomedicals.com **Product number: K8-3335**

Product name: Seta-555-NHS

General Data

Molecular Mass: 986.08 (protonated form)

Solubility: Water, Alcohol, DMF, DMSO **Insoluble:** Acetone, Chloroform, Toluene

Storage: Store in absence of light, desiccate and refrigerate

Description

• Seta-555-NHS (K8-3335) is highly hydrophilic, amine-reactive fluorescent label containing one reactive NHS-ester group with almost identical absorption and emission as Cy3™ or Alexa Fluor™ 555 and can therefore be used with these filter sets. It combines high photostability and brightness.

Applications

· Covalent labeling of proteins, amino-modified DNA and amino-modified oligonucleotides

Advantages


- · Perfectly suited for excitation with the 532-nm laser
- Good aqueous solubility: this label does not alter the solubility of the protein conjugate
- Low molecular weight: Seta-555 does not add substantial mass to the conjugates
- Photostability: Higher photostability as compared to Alexa Fluor™ 555 or Cy3
- Ideal for non-radioactive labeling of proteins, amino-modified DNA probes and amino-modified oligonucleotides

Spectral Data


Solvent System: phosphate buffer pH 7.4

Sample	Dye-to-protein Ratio	Absorption max. [nm]	Extinction Coefficient [M ⁻¹ cm ⁻¹]	Fluorescence max. [nm]	Quantum Yield ¹ [%]
Free dye	_	555	155,000	570	7
IgG conjugate 1	1.0	556		571	13
IgG conjugate 2	9.5	556		572	10


¹ Excitation at 490 nm



Absorption and emission spectrum of a **Seta-555** in phosphate buffer (pH 7.4)




Absorption and emission spectrum of a **Seta-555** — **IgG conjugate** in phosphate buffer (pH 7.4, Dye-to-protein ratio ~1.0)




http://www.setabiomedicals.com e-mail: info@setabiomedicals.com Product number: K8-3335

Product name: Seta-555-NHS



Quantum yield vs. dye-to-protein ratio (D/P) of **Seta-555** — **IgG** conjugates in phosphate buffer (pH 7.4) as compared to **Cy3** conjugates



Total brightness (Q.Y. $x \in x$ D/P) of **Seta-555** — **IgG conjugates** in phosphate buffer (pH 7.4) as compared to **Cy3** conjugates